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FOLEY CATHETERS - PERMANENT

PRODUCT DESCRIPTION

Catheters are designated for routine measurement of urinary bladder emptying, for drainage after surgery or urinary bladder irrigation.

PREPARATION

Take the catheter out from the package under septic conditions. Visual control is focused on potential damage to the catheter. In case of filling the balloon due to testing, it is necessary to empty it again before it is introduced. Use only sterile water.

INSTRUCTIONS FOR USE

Using lubricating substance or water-based gel, lubricate the catheter so that it is slippery (with one or more local anaesthetics as needed). Introduce the catheter into urethra under aseptic conditions. Make sure that the balloon lies behind the urinary bladder neck.

INSTRUCTIONS MANUAL

You can fill the balloon by a syringe of Luer type. Use sterile water. The balloon capacity is stated on the funnel. You must not overfill the balloon. Check the catheter throughput according to the back flow of urine in the drainage funnel. To unblock the balloon before removing the catheter place the syringe of Luer type on the valve and suck the balloon content slowly. You have to avoid excessive aspiration which can cause vacuum collapse of the filling lumen which would result in drainage failure.

Maximum period recommended for leaving the catheter in the bladder is 30 days. The inflatable balloon holds the catheter in the urinary bladder. If any fluid leaks from the balloon, e.g. soaking through the balloon membrane, then the catheter can fall out by mistake after some time. To prevent this you should suck out inflation fluid from the balloon every two weeks and refill it again according to the given nominal volume.

POSSIBLE COMPLICATION

If in exceptional cases the catheter cannot be unblocked, it is necessary to cut catheter pin with the valve on branching with sharp scissors. It is also possible to let the balloon burst. Make sure in this case that you remove the rest of it from the urinary bladder. Improperly placed catheter could cause damage to the urethra mainly if it is filled inside it.

WARNING

Use only for urological purposes. Use only water-based lubricants. The balloon is to be filled with sterile water or sterile water with glycerine in a ratio 9 : 1. Never use clamps - use closure (plugs) of urinary catheters. When taking the urinary sample, do not pierce the catheter with a needle; it could harm a patient or damage the catheter which would result in its replacement.

STORAGE

Do not store under extreme temperatures and in the damp environment.

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Keep away
from sunlight

Keep Dry

Caution, consult
accompanying documents

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